

DASYLab Serial Number Installation

In order to activate your DASYLab program then you will need to enter a valid serial number. The serial number will be provided after the purchase of the DASYLab program. The serial number can be used to activate either a 30-day trial installation or a new program installation. Use one of the following procedures below to activate your program.

Activate a 30-day Trial Installation

If you have previously installed a 30-day trial version of DASYLab and then purchased the program then you should use the following steps to enter your new serial number.

1. Click the Windows “Start” button in lower left corner of the Windows desktop.
2. Click on the “DASYLab 12.0” menu selection to expand and view the sub-menu selections.
3. Click on the “Configurator 12.0” sub-menu selection. See image 1.0 below.




Image 1.0

4. The Configurator 12.0 screen will open where the serial number can be entered. See image 2.0 below.
5. Ensure that the “Install Evaluation Version” checkbox is not checked.
6. Enter the serial number in the dialog box and press the “OK” button. The program may ask to reboot the computer in order to complete the process.
7. The process is now complete.

Activate a New Installation

During a new DASYLab program installation then you will be prompted to enter a valid program serial number. This serial number is supplied upon completion of a DASYLab program purchase. Use the following steps to enter your serial number.

1. Run the program installer, follow the prompts and answer the questions as they appear.
2. Upon arriving at the Configurator screen (see Image 2.0 below) enter the serial number and ensure that the “Install Evaluation Version” checkbox is not checked.
3. Enter your User Name and Company Name on the Configurator screen and press OK when done to continue with the installation.
4. The process will be complete when you have finished the remainder of the installation process.




Image 2.0

Contact Information:

DGH Corporation, USA.
P.O. Box 5638
Manchester, NH 03108
(603) 622-0452
www.dghcorp.com
information@dghcorp.com